

RELEASE MEMO FOR SOLABS QM 10.10.0-9

Status: Approved & Effective

Effective Date: 28-Jan-2026

Control Number: MEMO000074

Version: 1.0

DOCUMENT INFORMATION TABLE

Name	Release Memo for SOLABS QM 10.10.0-9
Document Type	Internal Memo
Description	This is a release memo that confirms that the corrective release 10.10.0-9 is qualified for use.
Control Number	MEMO000074
Version (Internal Version)	1.0 (1.0)
Effective Date	28-Jan-2026
Next Review Date	N/A
Additional Information	< empty >
Area	Product and Quality Management
ISO Reference	8. Operation
Originator (Author)	mboire@solabs.com (Martine Boire)
Document Owner	N/A
Document Coordinator	N/A
Unique ID	b02037b9-25ed-44d4-8cde-1d857c27fc9c

RELATED ITEMS

N/A

RELEASE MEMO FOR SOLABS QM 10.10.0-9

Status: Approved & Effective

Effective Date: 28-Jan-2026

Control Number: MEMO000074

Version: 1.0

APPROVAL TABLE

Status ¹	Selected for Approval username (Full Name), Title	Approved By username (Full Name)	Meaning	Approval Date
APP	mboire@solabs.com (Martine Boire), <i>Head of Product & Quality Management</i>	mboire@solabs.com (Martine Boire)	Originator (Author)	28-Jan-2026 06:06:54 UTC-05:00
APP	pgaudreau@solabs.com (Philippe Gaudreau), <i>CEO</i>	pgaudreau@solabs.com (Philippe Gaudreau)	Approver	28-Jan-2026 09:45:23 UTC-05:00

¹ Legend:

NS: Not Started

AWA: Awaiting Approval

APP: Approved

REJ: Rejected

STO: Stopped

From: Martine Boire, Head of Product and Quality Management

Date: January 27, 2026

Re: Release Memo for SOLABS QM 10.10.0-9 (Corrective Release)

Following the qualification of SOLABS QM10 release 10.10.0, an issue was identified related to System Attributes when an attribute name is modified:

System Attributes are created and linked to a document type. When documents of this type are created or modified, the associated System Attribute is displayed in the document metadata. However, if the name of a System Attribute is later changed, a display issue occurs. Documents created before the name change, both current and previous versions, no longer display the attribute in the document metadata, including the Document view page, Authoring History, and Audit Trail.

To address this issue, **corrective release 10.10.0-9** was issued and successfully qualified using the following test:

- **Test SOLABS_QMX_2024**

Testing Results:

- Test **SOLABS_QMX_2024** was executed successfully.
- No deviations were encountered during testing of corrective release **10.10.0-9**.

The successful execution of this test confirms that the issue has been resolved.

This memo confirms that **SOLABS QM 10.10.0-9** is qualified for deployment in a production environment and can be installed wherever **SOLABS QM 10.10.0** (or any other 10.10.0 corrective release) is in use.