

# RELEASE NOTES FOR SOLABS QM 10.11.0 (DRAFT – ISSUED DURING QUALIFICATION)

Internal Tracking (Ticket No.)	Area(s) Impacted	WHAT is New?	WHY?
CATEGORY: <u>UI/UX</u>			
Refer to list at the end of table*	All sections	<p>The user interface refresh introduced in Release 10.11 is an important step in the ongoing evolution of the SOLABS QM10 platform. As customer needs, technology standards, and user expectations continue to evolve, it became essential to modernize the user experience to ensure the system remains intuitive, efficient, and scalable. The updated interface is designed to improve usability, streamline navigation, and create a more consistent experience across the application.</p> <p>In addition to enhancing the overall user experience, the new UI provides a stronger foundation for future innovation. The modernized framework enables SOLABS to introduce new features more efficiently, maintain consistency across modules, and support continued growth of the platform while meeting current industry expectations for software usability and accessibility.</p> <p>Ultimately, the UI revamp reflects SOLABS' commitment to delivering a modern, user-centered quality management solution that helps customers work more effectively and confidently.</p>	
Refer to list at the end of table**	Search Function	<p>Complete review of Search Function</p> <p>Allow for wider adoption of the Search Function by providing filtering and ordering options, facilitating usage of supported search syntax, automatically saving search queries and including a <i>back to search results</i> option when leaving the search page, and narrowing fields "searched into" by default.</p>	
QMX-1545	Document	Document Status displayed in the document treeview and in the table views have a specific color per status.	User experience enhancement.
QMX-1701	Top Panel	A list of the 10 recently accessed folders is always available on the Top Panel.	Respond to client requests to facilitate direct and quick access to the document treeview.

## RELEASE NOTES FOR SOLABS QM 10.11.0 (DRAFT – ISSUED DURING QUALIFICATION)

Internal Tracking (Ticket No.)	Area(s) Impacted	WHAT is New?	WHY?
QMX-1707	Copy Link function	The Get Link menu is replaced by a Copy Link menu that automatically copies the link to your clipboard.	Save time, enhance user experience.
QMX-1721	QM APPs	A new process navigator is available in the process view page to navigate from one process to another (next, previous, first, last, next tenth, previous tenth).	Facilitate review and access to multiple quality processes during a single session.
QMX-1737	Top Panel	A link is displayed to the Root instance for end-users who have access to it. The Root instance is SOLABS' internal instance, which is used by system owners, administrators to approve documents, support periodic audits, etc.	Facilitates retrieval of the URL.
QMX-1738	Table Views	The possibility of going from one page to the other on the table views or expand/collapse the table is available at the top and the bottom of the tables.	Direct response to client requests.
QMX-1740	Views	Views have their specific URL, and Open Link in a new Tab feature is implemented for top menus, views, treeview.	Improved support to use multiple tabs within a single browser window.
QMX-1879	All Pages	An option to go back to the top of the page is available, the Copyright and Release information are displayed, and the client logo is displayed at the bottom of every system page.	Improve navigation.

## RELEASE NOTES FOR SOLABS QM 10.11.0 (DRAFT – ISSUED DURING QUALIFICATION)

Internal Tracking (Ticket No.)	Area(s) Impacted	WHAT is New?	WHY?
<b>QMX-2013</b>	Timeout Popup	The timeout popup is visible to users no matter where they are in the application i.e. if they are at the top or at the complete bottom of a form.	Direct response to client requests.
<b>QMX-2119</b>	Link Document Component	It is possible to load more than 100 results when linking a SOLABS' Document to another system item.	Direct response to client requests.
<b>QMX-2153</b>	Training	The name of the menu 'Training Status' is modified to 'My Training Activities'.	Direct response to client requests.
<b>CATEGORY: <u>QM APPS</u></b>			
<b>SQP-6877</b>	QM APP	P0007 – New Look and Feel	Harmonization with 10.11 look and feel.
<b>SQP-6878</b>	QM APP	P0044 – New Look and Feel	
<b>SQP-6928</b>	QM APP	P0052 – New Look and Feel	
<b>SQP-6934</b>	QM APP	P0036 – New Look and Feel	
<b>SQP-6940</b>	QM APP	P0050 – New Look and Feel	
<b>SQP-6941</b>	QM APP	P0055 – New Look and Feel	
<b>SQP-6956</b>	QM APP	P0034 – New Look and Feel	
<b>SQP-6957</b>	QM APP	P0035 – New Look and Feel	

## RELEASE NOTES FOR SOLABS QM 10.11.0 (DRAFT – ISSUED DURING QUALIFICATION)

Internal Tracking (Ticket No.)	Area(s) Impacted	WHAT is New?	WHY?
<b>SQP-6958</b>	QM APP	P0038 – New Look and Feel	
<b>SQP-6959</b>	QM APP	P0042 – New Look and Feel	
<b>SQP-6960</b>	QM APP	P0046 – New Look and Feel	
<b>SQP-6961</b>	QM APP	P0047 – New Look and Feel	
<b>SQP-6984</b>	QM APP	P0031 – New Look and Feel	
<b>SQP-6985</b>	QM APP	P0045 – New Look and Feel	
<b>SQP-6986</b>	QM APP	P0066 – New Look and Feel	
<b>CATEGORY: <u>CORRECTION OF KNOWN ISSUES</u></b>			
<b>QMX-572</b>	Maximum Character Validation	Method of Evaluation (Training Activity) and Comments (Task) fields display the maximum characters validation message.	Bug fix
<b>QMX-1725</b>	Training Notification	The Notification for Training Failed is correctly sent after the defined threshold that is now configurable. Default is set to 3.	Bug fix
<b>QMX-1732</b>	Training Confirmation with Assessment	Training can be recorded without errors when the activity is linked to an assessment.	Bug fix

## RELEASE NOTES FOR SOLABS QM 10.11.0 (DRAFT – ISSUED DURING QUALIFICATION)

Internal Tracking (Ticket No.)	Area(s) Impacted	WHAT is New?	WHY?
<b>QMX-1747</b>	Login	The awesome font message is never displayed when login out and login back in the application.	Bug fix
<b>QMX-1901</b>	Document - Collaborative Review	If there is a communication issue between QM10 and SharePoint while collaborating in the collaborative review mode, an error message is displayed to the user, and the system retries to communicate and access the document.	Bug fix
<b>QMX-1932</b>	Related Items - Expand / Collapse All	Expand All/Collapse All components are working correctly per Related Items table.	Bug fix
<b>QMX-1951</b>	Setup - Table Views	It is possible to access pages after the first page on the Document Workflow Templates table in the Setup Section.	Bug fix
<b>QMX-1988</b>	Login	When login in QM10 using a SOLABS' link, the user is directed to the concerned page and not the home page when using Azure AD and only one login is required.	Bug fix
<b>QMX-2022</b>	Link Document (Process, Task, Training)	Supporting Documents are filtered out when linking a SOLABS document to a Process, a Task or a Training.	Bug fix
<b>QMX-2069</b>	Secondary Tasks - QM Apps Waiting Phases	Process moves out of waiting when there are secondary tasks link using the new secondary tasks behaviors (restriction of privileges on secondary tasks)	Bug fix

## RELEASE NOTES FOR SOLABS QM 10.11.0 (DRAFT – ISSUED DURING QUALIFICATION)

Internal Tracking (Ticket No.)	Area(s) Impacted	WHAT is New?	WHY?
<b>QMX-2155</b>	Process - Modify Values	It is possible to use the Modify Values option when the process is completed or in a Waiting Phase.	Bug fix
<b>QMX-2196</b>	Document - Collaborative Review	It is possible to use the collaborative review function for users having capital letters in their User Principal Name in Azure directory.	Bug fix
<b>QMX-2197</b>	Calendars (all sections)	When a date is selected with a calendar between March 10 and November 02, the correct date with UTC -4:00 is displayed in the date field.	Bug fix
<b>QMX-2198</b>	Training - Create Curriculum	The list of documents when searching for a document to be linked as the curriculum document is loading correctly.	Bug fix
<b>CATEGORY: FEATURE-LOGIG CHANGE</b>			
<b>QMX-1880</b>	Setup - User Account	User account options for a user for his/her own account is available in the top Navigation bar instead of from the Setup section.	UI Enhancement.
<b>QMX-1924</b>	PDF Conversion	Adlib Job Ticket is accessed remotely and can be modified without impacting on the SOLABS QM10 Service.	Support enhancements to the Adlib Job Ticket without impacting clients.
<b>QMX-1931</b>	QM APPs - Secondary Tasks	Possibility of an assignee to act on Secondary Tasks when the process is completed.	

## RELEASE NOTES FOR SOLABS QM 10.11.0 (DRAFT – ISSUED DURING QUALIFICATION)

Internal Tracking (Ticket No.)	Area(s) Impacted	WHAT is New?	WHY?
<b>QMX-2049</b>	QM APPS - Secondary Tasks	Privileges set on secondary tasks related to any process in QM10 for users having either Act On or Act On & Reassign on process are now set to View.	Refer to details here: <a href="https://docs.solabs.com/hc/en-us/articles/35000066114455-Secondary-Tasks-Updated-Privilege-Mapping-Access-Restrictions-and-Configurable-Privilege-Matrix">https://docs.solabs.com/hc/en-us/articles/35000066114455-Secondary-Tasks-Updated-Privilege-Mapping-Access-Restrictions-and-Configurable-Privilege-Matrix</a>
<b>QMX-2050</b>	QM APPS - Secondary Tasks	Creating, modifying and deleting secondary tasks from within a QM APP is aligned with new revisit privileges on tasks.	
<b>QMX-2051</b>	QM APPS - Secondary Tasks	Secondary Tasks are displayed in the Related Items no matter which privileges the user has while In process and Completed. The Act On and the View menus are not available until the tasks are released.	
<b>QMX-2052</b>	QM APPS – Secondary Tasks	Any process step releasing secondary task can be confirmed by users no matter which privileges they have on the process.	
<b>QMX-2200</b>	Document - System Attribute	System attributes can be renamed without any impact on the attributes displayed in each document version.	
<b>CATEGORY: <u>SECURITY AND PERFORMANCE</u></b>			
<b>QMX-174</b>	QM APPS	Performance issues when many users are starting processes at the same time. This is fixed.	Performance improvement
<b>QMX-1592</b>	Database	A view is added in the database to get a human readable log of the system timers.	Performance improvement

## RELEASE NOTES FOR SOLABS QM 10.11.0 (DRAFT – ISSUED DURING QUALIFICATION)

Internal Tracking (Ticket No.)	Area(s) Impacted	WHAT is New?	WHY?
<b>QMX-1711</b> <b>QMX-2074</b> <b>QMX-2076</b> <b>QMX-2080</b> <b>QMX-2081</b> <b>QMX-2095</b> <b>QMX-2098</b> <b>QMX-2104</b>	Security	Security Vulnerability Remediations.	N/A

**NOTE:** The list above is not the full list of changes for release 10.11.0. Minor and technical changes (such as packaging or infrastructure enhancements) are not appearing in the release notes. Refer to the **FORM000011 – QM10.11.0 - Change List** for the complete list of changes.

**\*Internal Tracking (Ticket #):** QMX-1693, QMX-1694, QMX-1705, QMX-1745, QMX-1764, QMX-1775, QMX-1776, QMX-1804, QMX-1824, QMX-1825, QMX-1826, QMX-1833, QMX-1834, QMX-1835, QMX-1841, QMX-1849, QMX-1850, QMX-1851, QMX-1852, QMX-1853, QMX-1854, QMX-1855, QMX-1856, QMX-1857, QMX-1858, QMX-1859, QMX-1860, QMX-1861, QMX-1862, QMX-1863, QMX-1864, QMX-1865, QMX-1866, QMX-1867, QMX-1868, QMX-1869, QMX-1870, QMX-1871, QMX-1872, QMX-1873, QMX-1874, QMX-1875, QMX-1876, QMX-1877, QMX-1878, QMX-1879, QMX-1881, QMX-1882, QMX-1883, QMX-1884, QMX-1893, QMX-1894, QMX-1925, QMX-1939, QMX-1940, QMX-1941, QMX-1966, QMX-1968, QMX-1999, QMX-2004, QMX-2005, QMX-2014, QMX-2015, QMX-2016, QMX-2017, QMX-2035, QMX-2120, QMX-2121, QMX-2134, QMX-2141, QMX-2174

**\*\*Internal Tracking (Ticket #):** QMX-1716, QMX-1806, QMX-1808, QMX-1809, QMX-1809, QMX-1810, QMX-2011, QMX-2124, QMX-2170, QMX-2171, QMX-2173, QMX-2179.