

ONELOGIN/SOLABS QM - CONFIGURATION PARAMETERS FORM

Please provide a value for the following parameters and return SOLABS*		
Parameter	Section in administrative interface	Value
OneLogin SOLABS QM App Name	Portal -> Display Name	
OneLogin instance SAML URL	Enable SAML 2.0 -> SAML 2.0 Endpoint	
Region of the OneLogin instance	N/A	US for USA customers EU for European customers
SOLABS QM domain from the SOLABS app URL as configured in OneLogin	Application data -> ACS	
Client ID	Developer -> API Credentials	
Client Secret	Developer -> API Credentials	

*** Please make sure to also send the generated certificate file along with this form.**