

RELEASE MEMO FOR SOLABS QM 10.9.1-4

Status: Approved & Effective

Effective Date: 2024-02-01

Control Number: MEMO000058

Version: 1.0

DOCUMENT INFORMATION TABLE

Name	Release Memo for SOLABS QM 10.9.1-4
Document Type	Internal Memo
Description	This is a release memo that confirms that the corrective release 10.9.1-4 is acceptable for use.
Control Number	MEMO000058
Version (Internal Version)	1.0 (1.0)
Effective Date	2024-02-01
Next Review Date	N/A
Additional Information	< empty >
Area	Product and Quality Management
ISO Reference	8. Operation
Originator (Author)	mboire@solabs.com (Martine Boire)
Document Owner	N/A
Document Coordinator	N/A
Unique ID	ff719414-6bdb-4c83-8dad-92de1be9b85c

RELATED ITEMS

N/A

RELEASE MEMO FOR SOLABS QM 10.9.1-4

Status: Approved & Effective	Effective Date: 2024-02-01
Control Number: MEMO000058	Version: 1.0

APPROVAL TABLE

Status ¹	Selected for Approval username (Full Name), Title	Approved By username (Full Name)	Meaning	Approval Date
APP	kevin.b.langlois@solabs.com (Kevin Byrne-Langlois), Head of Engineering	kevin.b.langlois@solabs.com (Kevin Byrne-Langlois)	Approver	2024-02-01 15:42:36 UTC-05:00
APP	mboire@solabs.com (Martine Boire), Head of Product & Quality Management	pgaudreau@solabs.com * (Philippe Gaudreau)	Approver	2024-02-01 16:04:50 UTC-05:00

* Approved on behalf

¹ Legend:
NS: Not Started
AWA: Awaiting Approval
APP: Approved
REJ: Rejected
STO: Stopped

FORM000011 - SOLABS QM CORE SOFTWARE - CHANGE LIST 10.9.1-4

Status: Approved & Effective

Effective Date: 2024-02-01

Control Number: FORM000011 - QM10.9.1-4

Version: 1.0

DOCUMENT INFORMATION TABLE

Name	FORM000011 - SOLABS QM CORE SOFTWARE - CHANGE LIST 10.9.1-4
Document Type	Change List
Description	List of changes for Corrective Release 10.9.1-4.
Control Number	FORM000011 - QM10.9.1-4
Version (Internal Version)	1.0 (1.0)
Effective Date	2024-02-01
Next Review Date	N/A
Additional Information	< empty >
Originator (Author)	mboire@solabs.com (Martine Boire)
Document Owner	N/A
Document Coordinator	N/A
Unique ID	bbd83dbd-9448-475d-971b-b2964df5bdf8

RELATED ITEMS

N/A

FORM000011 - SOLABS QM CORE SOFTWARE - CHANGE LIST 10.9.1-4

Status: Approved & Effective

Effective Date: 2024-02-01

Control Number: FORM000011 - QM10.9.1-4

Version: 1.0

APPROVAL TABLE

Status ¹	Selected for Approval username (Full Name), Title	Approved By username (Full Name)	Meaning	Approval Date
APP	Any, Head of Engineering	kevin.b.langlois@solabs.com (Kevin Byrne-Langlois)	Approver	2024-02-01 15:49:18 UTC-05:00
APP	Any, Core Software Lead Developer	pdemers@solabs.com (Pascal Demers)	Approver	2024-02-01 15:28:09 UTC-05:00
APP	Any, Head of Product & Quality Management	schapdelaine@solabs.com * (Sophie Chapdelaine-Trépanier)	Approver	2024-02-01 15:53:46 UTC-05:00
APP	Any, Head of Client Success	emoore@solabs.com (Ericka Moore)	Approver	2024-02-01 15:53:50 UTC-05:00

* Approved on behalf

¹ Legend:

NS: Not Started

AWA: Awaiting Approval

APP: Approved

REJ: Rejected

STO: Stopped

SOLABS QM CORE SOFTWARE - CHANGE LIST – 10.9.1-4 (Corrective Release)



CONTROL NUMBER: FORM000011 – QM10.9.1-4

DOCUMENT VERSION: 1.0

DOCUMENT OWNER (Company): SOLABS

DOCUMENT OWNER (Department): Product and Quality Management

DOCUMENT AUTHOR(S): Martine Boire

Document History

Version	Date	Authors	Changes
1.0	Refer to document cover page	Martine Boire	First issue.

EFFECTIVE

TABLE OF CONTENTS

1 CHANGE LIST FOR CORE SOFTWARE SYSTEM RELEASE (Corrective Release)4

1.1 Change List for Release 10.9.1 to Release 10.9.1-44

EFFECTIVE

1 CHANGE LIST FOR CORE SOFTWARE SYSTEM RELEASE (CORRECTIVE RELEASE)

1.1 CHANGE LIST FOR RELEASE 10.9.1 TO RELEASE 10.9.1-4

Internal Tracking (Ticket No.)	Dev # (if applicable)	Description of Issue / Change	Test No.	Test Name	Status
CATEGORY: SECURITY AND PERFORMANCE ENHANCEMENT					
QMX-1751	N/A	Azure SSO Connector: MS Graph API calls to update the password in SOLABS active directory can be turned-off by configuration. And when enabled, such calls must log an exception when failing and become non-blocking.	SOLABS_QMX_1751	Verification of Login and E-signature with Azure AD Login Type.	Completed and Closed

From: Martine Boire, Head of Product and Quality Management

Date: February 01, 2024

Re: Release Memo for SOLABS QM 10.9.1-4 (Corrective Release)

Following the qualification of SOLABS QM10 release 10.9.1:

An issue was found in the Azure AD Connector (MS Graph API) causing intermittent connection issues, preventing the clients from being able to log into SOLABS QM10 or to provide an e-signature, sporadically. This issue applies only to clients using the Azure AD SSO Login Type.

A fix has been done to allow users to log into SOLABS QM10 even if there is an error on the MS Graph API side.

The corrective release 10.9.1-4 has then been issued for the qualification of the fix using the following test:

- Test SOLABS_QMX_1751

Testing Results:

- Test SOLABS_QMX_1751 was executed successfully.

There was no deviation generated during the testing of the corrective release 10.9.1-4.

The test has been executed with a successful result which confirms that the issue has been successfully corrected and implemented.

This memo confirms that SOLABS **QM 10.9.1-4** corrective release is qualified for use in a production instance and can be installed wherever the SOLABS QM 10.9.1 (or 10.9.1-1, 10.9.1-2 and 10.9.1-3) is installed.