

ONELOGIN/SOLABS QM - CONFIGURATION PARAMETERS FORM

Please provide a value for the following parameters and return SOLABS*		
Parameter	Section in administrative interface	Value
OneLogin SOLABS QM App Name	Info->Portal -> Display Name	
OneLogin instance SAML URL	SSO->Enable SAML 2.0 -> SAML 2.0 Endpoint	
Region of the OneLogin instance	N/A	US for USA customers EU for European customers
SOLABS QM10 app URL as configured in OneLogin	Configuration->Application details->ACS (Consumer) URL Validator*	
Reply URL as configured in OneLogin	Configuration->Application details->ACS (Consumer) URL*	
Client ID	Developer -> API Credentials	
Client Secret	Developer -> API Credentials	

* Please make sure to also send the generated certificate file along with this form.