

RELEASE MEMO FOR SOLABS QM 10.10.0-5 (CORRECTIVE RELEASE)

Status: Approved & Effective

Effective Date: 12-Jun-2025

Control Number: MEMO000066

Version: 1.0

DOCUMENT INFORMATION TABLE

Name	Release Memo for SOLABS QM 10.10.0-5 (Corrective Release)
Document Type	Internal Memo
Description	This is a release memo that confirms that the corrective release 10.10.0-5 is qualified for use.
Control Number	MEMO000066
Version (Internal Version)	1.0 (1.0)
Effective Date	12-Jun-2025
Next Review Date	N/A
Additional Information	< empty >
Area	Product and Quality Management
ISO Reference	8. Operation
Originator (Author)	mboire@solabs.com (Martine Boire)
Document Owner	N/A
Document Coordinator	N/A
Unique ID	4b0d5c99-8d24-4372-aa2b-15546d74a2f9

RELATED ITEMS

N/A

RELEASE MEMO FOR SOLABS QM 10.10.0-5 (CORRECTIVE RELEASE)

Status: Approved & Effective

Effective Date: 12-Jun-2025

Control Number: MEMO000066

Version: 1.0

APPROVAL TABLE

Status ¹	Selected for Approval username (Full Name), Title	Approved By username (Full Name)	Meaning	Approval Date
APP	mboire@solabs.com (Martine Boire), <i>Head of Product & Quality Management</i>	mboire@solabs.com (Martine Boire)	Originator (Author)	12-Jun-2025 14:38:13 UTC-04:00
APP	pgaudreau@solabs.com (Philippe Gaudreau), <i>CEO</i>	pgaudreau@solabs.com (Philippe Gaudreau)	Approver	12-Jun-2025 17:07:18 UTC-04:00

¹ Legend:

NS: Not Started

AWA: Awaiting Approval

APP: Approved

REJ: Rejected

STO: Stopped

From: Martine Boire, Head of Product and Quality Management

Date: June 12, 2025

Re: Release Memo for SOLABS QM **10.10.0-5** (Corrective Release)

Following the qualification of SOLABS QM10 release 10.10.0:

An issue was identified with the **DS_Document Metadata** report: when the number of System Attributes configured in the setup section exceeds a certain threshold (approximately 200), the report fails to execute. As a result, the **DS_Document_Metadata** report and any other reports relying on the same dataset cannot be generated, and attempting to access them results in an error. This issue was resolved in corrective release **10.10.0-1**.

Following the qualification of the SOLABS QM10 corrective release **10.10.0-1**, it was discovered that the deployment package inadvertently included Related Items tables that were not intended to be deployed. To address this, corrective release **10.10.0-5** was created to repackage the deployment without the unintended components and effectively removing the Related Items tables that had been deployed with version 10.10.0-1.

The **corrective release 10.10.0-5** was then issued and successfully qualified using the following tests:

- Test SOLABS_QMX_2006_1
- Test SOLABS_QMX_2006_2

Testing Results:

- Test SOLABS_QMX_2006_1 and SOLABS_QMX_2006_2 were executed successfully.
- No deviations were encountered during testing of corrective release **10.10.0-5**.

The successful execution of this test confirms that the issue has been resolved.

This memo confirms that **SOLABS QM 10.10.0-5** is qualified for deployment in a production environment and can be installed wherever **SOLABS QM 10.10.0** and **10.10.0-1** are in use.