

RELEASE MEMO FOR THE SOLABS QM 10.10.0-7

Status: Approved & Effective

Effective Date: 02-Oct-2025

Control Number: MEMO000068

Version: 1.0

DOCUMENT INFORMATION TABLE

Name	Release Memo for the SOLABS QM 10.10.0-7
Document Type	Internal Memo
Description	This is a release memo that confirms that the corrective release 10.10.0-7 is qualified for use.
Control Number	MEMO000068
Version (Internal Version)	1.0 (1.0)
Effective Date	02-Oct-2025
Next Review Date	N/A
Additional Information	< empty >
Area	Product and Quality Management
ISO Reference	8. Operation
Originator (Author)	mboire@solabs.com (Martine Boire)
Document Owner	N/A
Document Coordinator	N/A
Unique ID	453210f9-20e1-4875-afea-2fa1ee17270b

RELATED ITEMS

N/A

RELEASE MEMO FOR THE SOLABS QM 10.10.0-7

Status: Approved & Effective

Effective Date: 02-Oct-2025

Control Number: MEMO000068

Version: 1.0

APPROVAL TABLE

Status ¹	Selected for Approval username (Full Name), Title	Approved By username (Full Name)	Meaning	Approval Date
APP	mboire@solabs.com (Martine Boire), <i>Head of Product & Quality Management</i>	mboire@solabs.com (Martine Boire)	Originator (Author)	02-Oct-2025 10:02:44 UTC-04:00
APP	pgaudreau@solabs.com (Philippe Gaudreau), <i>CEO</i>	pgaudreau@solabs.com (Philippe Gaudreau)	Approver	02-Oct-2025 15:21:46 UTC-04:00

¹ Legend:

NS: Not Started

AWA: Awaiting Approval

APP: Approved

REJ: Rejected

STO: Stopped

From: Martine Boire, Head of Product and Quality Management

Date: October 01, 2025

Re: Release Memo for SOLABS QM **10.10.0-7** (Corrective Release)

Following the qualification of SOLABS QM10 release 10.10.0:

An issue was identified in the Training section when creating a Training Curriculum for Document Revision. When the system contains a large number of Supporting Documents (i.e., documents attached to tasks, trainings, or processes) along with numerous Retired or Superseded documents, performance issues may occur in the search component used to link Curriculum Documents. As a result, the document list may take an extended time to load, or in some cases, may fail to load entirely. This issue has been resolved in corrective release **10.10.0-7**

The **corrective release 10.10.0-7** was issued and successfully qualified using the following tests:

- **Test SOLABS_QMX_1998**

Testing Results:

- Test **SOLABS_QMX_1998** was executed successfully.
- No deviations were encountered during testing of corrective release **10.10.0-7**.

The successful execution of this test confirms that the issue has been resolved.

This memo confirms that **SOLABS QM 10.10.0-7** is qualified for deployment in a production environment and can be installed wherever **SOLABS QM 10.10.0, 10.10.0-1, 10.10.0-2, 10.10.0-3, 10.10.0-5** or **10.10.0-6** are in use.