

RELEASE MEMO FOR SOLABS QM 10.8.1-2

Status: Approved & Effective

Effective Date: 2020-04-17

Control Number: MEMO000023

Version: 1.0

DOCUMENT INFORMATION TABLE

Name	Release Memo for SOLABS QM 10.8.1-2
Document Type	Internal Memo
Description	This is the release meme that confirms that the corrective release 10.8.1-2 is acceptable for use.
Control Number	MEMO000023
Version (Internal Version)	1.0 (1.0)
Effective Date	2020-04-17
Next Review Date	N/A
Additional Information	< empty >
Area	Product and Quality Management
ISO Reference	8. Operation
Originator (Author)	mboire@solabs.com (Martine Boire)
Document Owner	N/A
Document Coordinator	N/A
Unique ID	e8a1a1df-2ee7-4d45-89f5-306e2f2b26d7

RELATED ITEMS

N/A

RELEASE MEMO FOR SOLABS QM 10.8.1-2

Status: Approved & Effective

Effective Date: 2020-04-17

Control Number: MEMO000023

Version: 1.0

APPROVAL TABLE

Status ¹	Selected for Approval username (Full Name), Title	Approved By username (Full Name)	Meaning	Approval Date
APP	kevin.b.langlois@solabs.com (Kevin Byrne-Langlois), <i>Head of Engineering</i>	kevin.b.langlois@solabs.com (Kevin Byrne-Langlois)	Approver	2020-04-17 13:52:28 UTC-04:00
APP	mboire@solabs.com (Martine Boire), <i>Head of Product & Quality Management</i>	schapdelaine@solabs.com * (Sophie Chapdelaine-Trépanier)	Approver	2020-04-17 14:56:07 UTC-04:00

* Approved on behalf

¹ Legend:

NS: Not Started

AWA: Awaiting Approval

APP: Approved

REJ: Rejected

STO: Stopped

EFFECTIVE

To: File
From: Martine Boire, Head of product and Quality Management
Date: April 17, 2020
Re: Release Memo for SOLABS QM 10.8.1-2 (Corrective Release)

Following the qualification of SOLABS QM10 release 10.8.1, a modification to a new behavior has been asked in regards to the new Timezone feature and the display on the document cover pages (UTC display). To this has been added a change in the system configurations to include an option to modify the document cover pages look and feel. Thus, a new look and feel for the cover pages is now available. Finally, a fix to the Assign by End User page (Training Administration section) where the calendars were not working anymore has been included.

These have been tested using the two (2) following tests scripts.

- For the testing of the modification related to the new timezone implemented in 10.8.1 release (UTC display in cover pages) and the new cover pages look and feel, the following test was executed:
 - o Test SOLABS_QMX_1237_QMX_1252.
- For the fix to the Assign by End User Page Calendars, the following test was executed:
 - o Test SOLABS_QMX_1263

Testing Results:

- o Test SOLABS_QMX_1237_QMX_1252 was executed successfully
- o Test SOLABS_QMX_1263 was executed successfully

There was no deviation generated during the testing of the corrective release 10.8.1-2.

This memo confirms that SOLABS QM 10.8.1-2 is acceptable for use in a production instance and can be installed wherever the SOLABS QM V10.8.1 is installed.

ANNEX 1: CHANGE LIST FOR CORRECTIVE RELEASE 10.8.1-2 (CL-CR-SOL-10.8.1-2)

CHANGE LIST FOR VERSION 10.8.1 TO VERSION 10.8.1-2

Internal Tracking (Ticket No.)	Dev # (if applicable)	Description of Issue / Change	Instruction for testing	Status
QMX-1237	N/A	Improvements to the Document Metadata page's look and feel.	Execute the following test: SOLABS_QMX_1237_QMX-1252	Formally tested, successful.
QMX-1252	N/A	Date/Timestamp issues in Document Metadata and Headers following introduction of time zones.	Execute the following test: SOLABS_QMX_1237_QMX-1252	Formally tested, successful.
QMX-1263	N/A	The Calendars in the Assign by End User page (Training Admin) do not work anymore	Execute the following test: SOLABS_QMX_1263	Formally tested, successful.