

SOLABS  
 777, de la Commune West  
 Suite 330  
 Montreal, Quebec  
 H3C 1Y1, Canada  
 514.287.1368 ext. 9  
 www.solabs.com



## RELEASE NOTES SOLABS QM 10.8.2

**Status:** Approved & Effective

**Effective Date:** 2020-11-30

**Control Number:** N/A

**Version:** 1.0

### DOCUMENT INFORMATION TABLE

|                                   |                                      |
|-----------------------------------|--------------------------------------|
| <b>Name</b>                       | Release Notes SOLABS QM 10.8.2       |
| <b>Document Type</b>              | Release Notes                        |
| <b>Description</b>                | < empty >                            |
| <b>Control Number</b>             | N/A                                  |
| <b>Version (Internal Version)</b> | 1.0 (1.0)                            |
| <b>Effective Date</b>             | 2020-11-30                           |
| <b>Next Review Date</b>           | N/A                                  |
| <b>Originator (Author)</b>        | gdixson@solabs.com (Gretchen Dixon)  |
| <b>Document Owner</b>             | FCT_Head of Client Success           |
| <b>Document Coordinator</b>       | N/A                                  |
| <b>Unique ID</b>                  | a1f741c9-ac47-4753-ae55-f3c7d8f3eb04 |

### RELATED ITEMS

N/A

SOLABS  
777, de la Commune West  
Suite 330  
Montreal, Quebec  
H3C 1Y1, Canada  
514.287.1368 ext. 9  
www.solabs.com



## RELEASE NOTES SOLABS QM 10.8.2

**Status:** Approved & Effective

**Effective Date:** 2020-11-30

**Control Number:** N/A

**Version:** 1.0

### APPROVAL TABLE

| Status <sup>1</sup> | Selected for Approval<br>username (Full Name), Title                                  | Approved By<br>username<br>(Full Name) | Meaning  | Approval Date                       |
|---------------------|---|--|----------|-------------------------------------|
| APP                 | pdemers@solabs.com (Pascal Demers),<br><i>Core Software Lead Developer</i>            | pdemers@solabs.com<br>(Pascal Demers)  | Approver | 2020-11-30<br>11:33:29<br>UTC-05:00 |
| APP                 | mboire@solabs.com (Martine Boire),<br><i>Head of Product &amp; Quality Management</i> | mboire@solabs.com<br>(Martine Boire)   | Approver | 2020-11-30<br>17:34:03<br>UTC-05:00 |

<sup>1</sup> Legend:

NS: Not Started

AWA: Awaiting Approval

APP: Approved

REJ: Rejected

STO: Stopped

EFFECTIVE

## RELEASE NOTES FOR SOLABS QM 10.8.2

| Internal Tracking (Ticket No.)                                     | Area(s) Impacted             | WHAT is new?  | WHY?   |
|--|------------------------------|---|--|
| <b>CATEGORY: <u>UI ENHANCEMENTS / ENHANCED USER EXPERIENCE</u></b> |                              |   |  |
| QMX-1317   | Document Cover Pages         | New look and feel available for Document Cover Pages.   | It is now possible to select between 2 different options for the look and feel of the Document Cover Pages on PDF Rendering Templates. The information remains the same; only the style is different.<br>This impacts only the document cover pages.<br><i>This is the integration in 10.8.2 of the change qualified in 10.8.1-2 as a corrective release for impacted clients.</i> |
| <b>CATEGORY: <u>NEW FEATURE</u></b>                                |                              |   |  |
| QMX-1275   | Login Page E-Signature       | Support for new Multi Factor Authentication (MFA) feature for clients using OneLogin SSO as the Login Type. | To allow use of multi factor authentication using a device such a cell phone or tablet with the use of Google Authenticator or OneLogin protect.   |
| QMX-1316   | Setup Login Page E-Signature | Integration with Microsoft Azure SSO as a Login Type option.  | To allow the integration of SOLABS with Microsoft Azure SSO.<br><i>This is the integration in 10.8.2 of a change qualified in 10.8.1-1 as a corrective release for impacted clients.</i>   |
| <b>CATEGORY: <u>CORRECTION OF KNOWN ISSUES</u></b>                 |                              |   |  |
| QMX-1318   | Document Cover Pages         | Remove the UTC stamp that was appended to the dates displayed in the header of Document Cover Pages.        | This information is not required for dates in the document header.<br><i>This is the integration in 10.8.2 of a change qualified in 10.8.1-2 as a corrective release for impacted clients.</i>   |

## RELEASE NOTES FOR SOLABS QM 10.8.2

| Internal Tracking (Ticket No.) | Area(s) Impacted | WHAT is new?   | WHY?   |
|--------------------------------|------------------|--|--|
| QMX-1319                       | Training Section | Calendar is now enabled to apply Due Dates in the Assign by End User option of the Training Administration menu.   | An issue was preventing use of the calendars in the Assign by End User option of the Training Administration menu. It was therefore not possible to apply Due Dates when assigning multiple Training Activities to a new user.<br><i>This is the integration in 10.8.2 of a change qualified in 10.8.1-2 as a corrective release for impacted clients.</i> |
| QMX-1320                       | Document Section | Correct an issue in corrective release 10.8.1-1 for Azure SSO clients using the Collaborative Review feature that was preventing the Review Cycle or Approval Cycle to be started. | The start of the Review and Approval Cycle function in systems that were updated with corrective release 10.8.1-1 was not working properly.<br><i>This is the integration in 10.8.2 of a change qualified in 10.8.1-4 as a corrective release for impacted clients.</i>  |
| QMX-1321                       | Training Section | Correct an issue with the Questionmark OnDemand integration to allow failing attempts to be properly recorded and reassigned.  | Failed assessments could not be reassigned for another try.<br><i>This is the integration in 10.8.2 of a change qualified in 10.8.1-4 as a corrective release for impacted clients.</i>  |
| QMX-1322                       |                  | Correct an issue in the SOLABS Active Directory so that passwords are correctly updated When using Azure SSO without SAML.   | When using Azure SSO without SAML, the password was not always correctly updated in the SOLABS Active Directory in some cases.<br><i>This is the integration in 10.8.2 of a change qualified in 10.8.1-4 as a corrective release for impacted clients.</i>   |