

## RELEASE NOTES SOLABS QM 10.9.1

**Status:** Approved & Effective

**Effective Date:** 2022-03-10

**Control Number:** N/A

**Version:** 2.0

### DOCUMENT INFORMATION TABLE

<b>Name</b>	Release Notes SOLABS QM 10.9.1
<b>Document Type</b>	Release Notes
<b>Description</b>	< empty >
<b>Control Number</b>	N/A
<b>Version (Internal Version)</b>	2.0 (2.0)
<b>Effective Date</b>	2022-03-10
<b>Next Review Date</b>	N/A
<b>Reason for Change</b>	Remove new macro-enabled Excel Report under QM APPs section since, per E. Moore, it will be released later than the updated SSRS Training Reports.
<b>Originator (Author)</b>	gdixson@solabs.com (Gretchen Dixon)
<b>Document Owner</b>	FCT_Head of Client Success
<b>Document Coordinator</b>	N/A
<b>Unique ID</b>	2cc2cc4f-be72-4788-aec1-b8013cf1daed

### RELATED ITEMS

N/A

SOLABS  
 9160 boul. Leduc  
 Suite 410  
 Brossard, Quebec  
 J4Y 0E3, Canada  
 1.877.322.1368  
 www.solabs.com



## RELEASE NOTES SOLABS QM 10.9.1

**Status:** Approved & Effective

**Effective Date:** 2022-03-10

**Control Number:** N/A

**Version:** 2.0

### APPROVAL TABLE

Status <sup>1</sup>	Selected for Approval username (Full Name), Title	Approved By username (Full Name)	Meaning	Approval Date
APP	mboire@solabs.com (Martine Boire), <i>Head of Product &amp; Quality Management</i>	mboire@solabs.com (Martine Boire)	Approver	2022-03-10 14:13:53 UTC-05:00
APP	<i>Any, Head of Client Success</i>	emoore@solabs.com (Ericka Moore)	Approver	2022-03-10 15:58:32 UTC-05:00

<sup>1</sup> Legend:

NS: Not Started

AWA: Awaiting Approval

APP: Approved

REJ: Rejected

STO: Stopped

EFFECTIVE

**RELEASE NOTES FOR SOLABS QM 10.9.1**

Internal Tracking (Ticket No.)	Area(s) Impacted	WHAT is new?	WHY?
<b>CATEGORY: <u>UI ENHANCEMENTS / ENHANCED USER EXPERIENCE</u></b>			
QMX-1285 QMX-1425 QMX-1426 QMX-1427 QMX-1428 QMX-1429	Reports	Update the following SSRS Training Reports to include additional parameters to filter by Assignment Type (Job Profile, End User or both), Department or No Department, to filter out assignments with future due dates. <ul style="list-style-type: none"> <li>• Employee Training Compliance</li> <li>• Training Compliance Distribution</li> <li>• Training Compliance User Distribution by Range</li> <li>• Training Compliance by Department</li> <li>• Departmental Training Compliance by User</li> <li>• Training Compliance Department Status by User</li> </ul>	Continuous improvement
<b>CATEGORY: <u>NEW FEATURE</u></b>			
N/A			
<b>CATEGORY: <u>QM APPS</u></b>			
N/A			
<b>CATEGORY: <u>CORRECTION OF KNOWN ISSUES</u></b>			
QMX-1448	Reports	Training Assignment as Trainer will be correctly displayed in the Employee Training History Report.	When Training Activity assignments were made to Trainers, they were not being displayed correctly.

**RELEASE NOTES FOR SOLABS QM 10.9.1**

Internal Tracking (Ticket No.)	Area(s) Impacted	WHAT is new?	WHY?
QMX-1459	Reports	Special character (&) will be displayed correctly in the Employee Training History Report.	The “&” character was not displaying properly for the Department names title over the table.
QMX-1461	Document	Multiple documents can be opened and closed one after the other from the Document Treeview.	After one document was opened and closed, another could not be opened from the Treeview without either first clicking on the View option or going out of the Treeview and back in again.
QMX-1472	Training	Questionmark OnDemand Assessments cannot be copied to a client folder unless it is from the Public folder or from another folder of the same client.	Must not be possible to link an assessment from folders other than Public or the logged in client folders.
QMX-1480	Document	Documents imported with past effective dates and periodic review requirements can be approved without issue.	An issue was preventing document approval if the Effective Date field was populated with a past date and there was a Periodic Review period set.
QMX-1482	Reports	The <i>isAssignedAs</i> column on the DS_Training_Users_Profiles dataset will be properly populated.	DS_Training_Users_Profiles was missing a value in the column <i>isAssignedAs</i> when the user was the Author or the Approver of the Curriculum Document.